

NASAL DRUG DELIVERY VIA THE ORAL ROUTE USING A PMDI

In this article, Laurent Vecellio, PhD, Research Engineer, University of Tours, and Scientific Director, Nemera (previously employed by Aerodrug-DTF Medical), D borah Le Pennec, Research Technician, University of Tours, and Alain Regard, Technology Product Manager, Nemera, discuss a study Nemera has funded into the retronose concept, using a pMDI to deliver to the nasal cavity via the oral route during exhalation through the nose.

INTRODUCTION

Nasal drug delivery is a non-invasive method that allows for a rapid, high and local therapeutic effect. It offers significant opportunities for new drug development looking to deliver systemic drugs, vaccines and treatments for the central nervous system. A recent study in patients with chronic rhinosinusitis has shown how deposition of corticosteroids in the nasal cavities can have an impact on clinical outcomes.¹ This study demonstrated the importance of the delivery device on drug efficacy.

Standard nasal sprays have limitations regarding the reproducibility and the deposition efficacy in the distal region of interest in the nasal cavities. An alternative, the "retronose" concept (Figure 1) has been proposed as a means to reduce variability and improve drug deposition in specific nasal zones, consisting of drug




Figure 1: The retronose concept.

"The "retronose" concept has been proposed as a means to reduce variability and improve drug deposition in specific nasal zones, consisting of drug administration through the buccal cavity during the nasal expiratory phase."

administration through the buccal cavity during the nasal expiratory phase.

Using this method, a different nebuliser concept has been developed for better drug deposition in the distal region of the nose,^{2,3} without lung deposition. Drug particles enter the nasal cavities through the rhinopharynx, which has a significant impact on drug deposition. Additionally, in a recent study, five asthmatics with rhinosinusitis were successfully treated with an aerosol therapy exhaled through the nose⁴ using a similar concept.

In the study funded by Nemera presented here, the use of a pressurised metered dose inhaler (pMDI), as an alternative to a nebuliser, for delivering drugs to the nose via the buccal cavity was explored. Specifically, the study focused on the influence of particle size and the expiratory flow rate on particle deposition in an upper airways model using a standard pMDI.



Dr Laurent Vecellio
Research Engineer
University of Tours, France
Scientific Director
Nemera
E: vecellio@med.univ-tours.fr



D borah Le Pennec
Research Technician
University of Tours, France



Alain Regard
Technology Product Manager
Nemera
E: alain.regard@nemera.net

Nemera
20, Avenue de la Gare - B.P. 30
38292 La Verpilli re Cedex
France

www.nemera.net




Figure 2: Experimental set up using the VCU upper airways model.

METHOD

A pMDI filled with HFA 134a propellant (no surfactant) was used with a 90 μ L valve system (Inhalia[®], Nemera) and an actuator (NM200, Bepak, Germany). Three different particle sizes (3, 12 and 20 μ m in terms of volume mean diameter) of a model drug were put in the canisters, resulting in three different pMDI suspensions (pMDI-A, pMDI-B, pMDI-C) delivering 100 μ g of drug per dose. Aerosol particle size produced was measured using a cascade impactor operating at 30 L/min (Next Generation Impactor, Copley Scientific, UK).

Aerosol deposition in the upper airways (Figure 2) was studied using an anatomical model⁵ developed by the Virginia Commonwealth University (Richmond, VA, US). The trachea model was connected to an absolute filter (L) and a humidified air source at three different flow rates: 30, 60 and 90 L/min. A vacuum pump connected to an absolute filter (E) was

located near to the nose model for collecting the totality of the exhaled aerosol from the model. Eight regions of interest were defined in the upper airway model:

- Mouth (M)
- Trachea (T)
- Oropharynx (O)
- Rhinopharynx (R)
- Upper part of the nasal cavity (U)
- Middle part of the nasal cavity (M)
- Bottom part of the nasal cavity (B)
- Nostrils (N).

The active compound was assayed by a spectrophotometric method.

RESULTS

The aerosol particle sizes produced by pMDI-A and pMDI-B, measured by cascade impaction, were characterised by a mass median aerodynamic diameter (MMAD) as shown in Table 1. The MMAD produced by pMDI-C

	pMDI-A	pMDI-B	pMDI-C
MMAD (μ m)	3.7 \pm 1.3	14.8 \pm 0.4	N/A
Particle size lower than 10 μ m (%)	70 \pm 17	29 \pm 8	3 \pm 1

Table 1: Aerosol particle size data for the three pMDIs.

“The concept of nasal drug delivery via the oral route using a pMDI in an upper airways model has been demonstrated *in vitro*, with promising results.”

could not be calculated due to its high deposition in the induction port of the cascade impactor. Regarding the percentage of particle size lower than 10 μ m. The results obtained with pMDI-C are consistent with the particle size put in the canister (Table 1).

An increase of drug penetration into the nasal cavity was observed with a decrease of particle size. No active compound was detected depositing into the filter corresponding to the lung model (L). No statistical difference was observed on the influence of expiratory flow rate on aerosol deposition in the oropharynx, mouth or trachea, and drug penetration into the nasal cavity ($p > 0.05$, Friedman test, GraphPad Prism V5) for all three flow rates were examined.

The drug deposited homogeneously in all regions of the nasal cavity (Figure 3). An exception was the upper part of the nasal cavity, where deposition was relatively low (less than 2%), independent of particle size or expiratory flow rate. An increase in particle size was associated with a decrease in exhaled fraction. An increase in expiratory air flow rate was also associated with a decrease in exhaled fraction. There was no statistical influence of expiratory flow rate on deposition fraction in the other regions of the nasal cavity, N, U, M, B and R ($p > 0.05$, Friedman test, GraphPad Prism V5).

CONCLUSION

The concept of nasal drug delivery via the oral route using a pMDI in an upper airways model has been demonstrated *in vitro*, with promising results. All anatomical regions, except for the upper part of the nasal cavity, were successfully targeted, with relatively homogenous deposition.

This nasal drug delivery system could be of interest for both local and systemic drug delivery, and for the delivery of vaccines.




Figure 3: Deposition in the nasal cavity model expressed in term of total drug delivered in the nasal cavity (N+U+M+B+E+R) for 3 μm (A), 12 μm (B) and 20 μm (C) (n=3).

ACKNOWLEDGEMENT

The authors would like to thank Dr P Worth Longest of Virginia Commonwealth University (Richmond, VA, US) for providing the upper airways model.

ABOUT THE COMPANY

Nemera designs, develops and manufactures nasal, buccal, auricular, ophthalmic, pulmonary, parenteral (passive safety devices autoinjectors, pens, and implanters), dermal

and transdermal drug delivery devices for the pharmaceutical, biotechnology and generics industries.

REFERENCES

1. Reychler G et al, "Effect of three-drug delivery modalities on olfactory function in chronic sinusitis". *Laryngoscope*, 2015, Vol 125(3), pp 549–555.
2. Kobayashi Y et al, "A novel therapeutic use of HFA-BDP metered dose inhaler

for asthmatic patients with rhinosinusitis: Case series". *Int J Clin Pharmacol Ther*, 2014, Vol 52(10), pp 914–919.

3. Sutton S, "Navigating the Nose". *The Medicine Maker*, Oct 2017.
4. Vecellio L, Chantrel G, Massardier JM, "Device for oral administration of an aerosol for the rhinopharynx, the nasal cavities or the paranasal sinuses". Patent WO2011080473.
5. Longest PW, Tian G, Hindle M, "Improving the Lung Delivery of Nasally Administered Aerosols During Noninvasive Ventilation—An Application of Enhanced Condensational Growth (ECG)". *J Aerosol Med Pulm Drug Deliv*, 2011, Vol 24(2), pp 103–118.

ABOUT THE AUTHORS


Laurent Vecellio is Scientific Director at Nemera, and research engineer in the Research Centre for Respiratory Diseases at the University of Tours, France. He obtained his PhD in 2002 and his accreditation to supervise university research in 2007. He worked for Aerodrug-DTF Medical in nebulisation from 2000 to 2018.

Déborah Le Pennec is a laboratory technician in the aerosol therapy team at the Research Centre for Respiratory Diseases, University of Tours, France. She works on different *in vitro* studies such as aerosol metrology, aerosol deposition, bioequivalence, and scintigraphy, and is experienced in biological and chemical analytical techniques.

Alain Regard, Technology Product Manager, Nemera, graduated with a degree in Polymer Engineering and Processing from ESP in Oyonnax, France. After a long experience in design and development in the automotive industry, he joined the company in 2010 as a product development leader. Mr Regard, today one of the key technical experts of Nemera's Innovation Center for Devices (ICD), leads the nasal and dermal developments. He drives some of Nemera's own IP projects as well as working on several customer product developments in the field of nasal and dermal applications.

Nemera

we put
patients
First



Advancia®

The advanced level
of patient adherence



Advancia® user-independent nasal pump guarantees a full dose delivery for every single spray. It is suitable for preserved and preservative-free formulations.

www.nemera.net • information@nemera.net • +33 (0)4 74 94 06 54

• INHALATION • EAR NOSE THROAT • PARENTERAL • OPHTHALMIC • DERMAL