



ACCELERATING TIME TO MARKET FOR NEW INJECTION DEVICES USING A NOVEL FORCE-MODELLING METHOD AND SOFTWARE

 **SANNER GROUP**
Protecting Health.



Bradley Sawyer, Elena Guss Tarazona and Matthew Latham, all at **Sanner Group's** Design Center of Excellence, **Springboard**, and **Charlie Bowen, Dr Jay Sayed, and Harriet Field**, all at **Pfizer**, describe new software co-developed by the companies to translate syringe test data into injection device performance, and to work backwards from device performance after shelf-life expiration to the maximum forces they can allow in a syringe test.

New injection devices are being asked to deliver larger volume injections, with higher viscosity formulations and increasing user expectations for device aesthetics, simplicity and quality. All injection devices risk performing inadequately when challenged with components, volumes and viscosities at the upper end of their tolerances. This risk of failure is present when developing devices for biosimilars and generics too.

As such, it is becoming ever more important to predict a device's behaviour early in development, as well as to predict the effect of shelf-life ageing. However, with increased drug viscosity and volume comes more complex force requirements, and basic methods of predicting plunger forces are not sufficiently capable. The nightmare scenario is that an injection device is designed, taken through tooling, automation and design verification testing, and then fails its performance requirements at the end of its shelf life.

Sanner Group's Design Centre of Excellence (Springboard) has developed new software that allows the translation of basic syringe test data into injection device performance. The software also enables engineers to work backwards from device performance after shelf life to

the maximum forces they can allow in a syringe test. This can be done early in development to avoid failures late in the programme.

WHY USE BREAK-LOOSE EXTRUSION FORCE DATA?

Break-loose and extrusion force (BLEF) tests are a common way to characterise prefilled syringes with plungers. In a BLEF test, a tensometer drives a plunger at constant speed and measures the (varying) force as the plunger moves down the syringe.

The force-distance chart created in a BLEF test provides useful information – the initial break-loose force of the stationary plunger, the extrusion force (which comprises the dynamic plunger-syringe friction and the drug formulation’s back-pressure) and a small dip in force indicative of an air bubble expelled after drug delivery (Figure 1).

ISO 11040-8 mandates that such tests use the complete, final system “as intended for use” and that the designer considers how forces vary with ageing and environmental factors. BLEF testing is valuable for characterising the relationships between forces and plunger speed, syringe or needle dimensions, and fluid properties.




Figure 1: The force-distance chart for a BLEF test.




Figure 2: Steps taken to generate BLEF predictive data.

However, it cannot be used to accurately predict injection times for a device because it does not reflect dynamic injection behaviour. Therefore, designers need a tool that can use common BLEF test data to provide accurate predictions for injection device performance.

To address this, Springboard has co-developed advanced modelling software

with Pfizer’s Devices Centre of Excellence. This software enables device developers to easily forecast injection times, anticipate variation due to ageing and set meaningful specifications early in the development process. The software can bridge the gap between constant-speed BLEF data and dynamic injection device performance without the need for guesswork, reducing the number of slow and costly prototype testing cycles and rework (Figure 2).

A NEW MODEL AND SOFTWARE FOR DEVICE PERFORMANCE

Software users simply upload their own BLEF data into the software and, in a few steps, can generate:

1. A predicted distribution of injection times for their specified device
2. Maximum extrusion force requirement predictions
3. A report for their design history file.

The model behind the software uses a time-stepped approach to break down an injection into discrete, calculable steps (Figure 3). Each loop begins with a known plunger position, which relates to the device’s spring length. From this, the spring force is known and is assumed to be approximately equal to the current extrusion force.




Figure 3: A time-stepped approach to break down an injection.

Using relationships from the user's BLEF data for the modelled syringe, the extrusion force is used to work back to the plunger speed and internal syringe pressure. By calculating the Hagen-Poiseuille equation across the needle, the pressure is related to the velocity of fluid ejected.

The plunger movement can then be calculated due to the reduction in injection volume, while also accounting for any changes in any air bubble volume due to pressure changes – this step is important to capture, but is complex and challenging to address through other predictive methods. The loop is repeated until the full volume has been delivered, with the total time elapsed giving the device injection time.

This software uses the Monte Carlo method to sample the random variation of all input parameters, which include:

1. User-specified tolerances on input parameters
2. The variance of real BLEF data.

COPING WITH DEVICE AGEING

Various factors lead to device performance degradation over time, most notably desiliconisation of the syringe, which leads to increased plunger-syringe friction. After shelf-life expiration, device performance is more likely to fail requirements such as injection time than on the day the device was manufactured. This is bad news for development projects because issues might not be discovered until years into a programme.

**"VARIOUS FACTORS
LEAD TO DEVICE
PERFORMANCE
DEGRADATION
OVER TIME,
MOST NOTABLY
DESILICONISATION
OF THE SYRINGE,
WHICH LEADS TO
INCREASED PLUNGER-
SYRINGE FRICTION."**

Therefore, Sanner Group and Pfizer used newly siliconised and fully desiliconised syringes as best and worst cases for plunger-syringe friction and used desiliconisation as a surrogate for ageing. In practice, while it is safe to say that syringes effectively move away from newly siliconised performance as they age, most will never get near to fully desiliconised, depending on their ageing properties.

The software has been extended to establish ageing trends using a small number of single-speed BLEF tests on aged syringes. This effect is modelled by increasing theoretical BLEFs until devices stall or breach injection time limits. This results in establishing the worst case

**"THE SOFTWARE
HAS BEEN EXTENDED
TO ESTABLISH
AGEING TRENDS."**

syringes that can be allowed. Then, by back-calculating from end-of-shelf life performance to initial conditions, the method can define a maximum allowable BLEF result at time zero. This allows engineers to specify the maximum forces allowed on a BLEF test at time zero where they can still be confident that the injection device would perform after shelf-life expiration.



Bradley Sawyer

Bradley Sawyer is a mechanical engineer at Sanner Group's Design Center of Excellence, Springboard, with particular experience in medical device concept design and injection device verification. He holds an MEng in mechanical engineering from the University of Cambridge (UK).

E: bradley.sawyer@springboard.pro



**Elena Guss
Tarazona**

Elena Guss Tarazona holds a BEng in mechanical engineering from the University of Warwick (UK) and an MRes in medical physics and biomedical engineering from University College London (UK). While at Sanner Group's Design Center of Excellence, Springboard, she has worked on a range of medical devices from early-stage R&D through to forensic engineering, with a focus on feasibility and rapid prototyping.

E: elena.gusstarazona@springboard.pro



**Matthew
Latham**

Matthew Latham, an engineer at Sanner Group's Design Center of Excellence, Springboard, has a background in mechanical and aerothermal engineering. He has worked on a broad range of medical devices across all stages of product lifecycles, with a particular focus on primary packaging, design verification and injection devices.

E: matthew.latham@springboard.pro

Springboard Ltd – a Sanner Group Company

Platinum Building, St John's Innovation Centre, Cowley Road,
Cambridge, CB4 0WS, United Kingdom
www.sanner-group.com

PRECISION ENGINEERED – SEAMLESSLY DELIVERED

As user expectations for device aesthetics, simplicity and quality continue to rise, the pressure on the medical device industry intensifies. Sanner empowers pharmaceutical and healthcare companies to meet these demands by uniting world-class design expertise, global manufacturing capabilities, regulatory know-how and

operational excellence. From next-generation injection systems and drug delivery devices to innovative diagnostics and connected health solutions, Sanner delivers a comprehensive, single-source approach that accelerates time to market. Built on a foundation of quality, trust and shared purpose, Sanner stands ready to shape the future of medical device development in close partnership with its customers.

SUMMARY

If injection time performance is not properly predicted (including predictions for after shelf life), a device programme risks substantial cost and time overruns. The cost and time come from having to implement time-consuming and expensive component modifications and reverification testing.

Instead, device developers can use the new model described here to change the physical build-and-test loop into a mathematical model, which specifies acceptable force limits for easy BLEF tests early in the programme. These predictive tools will become increasingly essential for efficient and reliable device development as drug delivery systems evolve to accommodate more demanding formulations. The modelling software offers a practical solution by translating BLEF test data into more accurate forecasts of injection performance over a device's lifespan.

“DEVICE DEVELOPERS CAN USE THE NEW MODEL DESCRIBED HERE TO CHANGE THE PHYSICAL BUILD-AND-TEST LOOP INTO A MATHEMATICAL MODEL, WHICH SPECIFIES ACCEPTABLE FORCE LIMITS FOR EASY BLEF TESTS EARLY IN THE PROGRAMME.”



Charlie Bowen

Charlie Bowen is a mechanical engineer at Pfizer, with a passion for developing medical devices. For the past three years, Mr Bowen has been developing combination products for Pfizer at the Devices Centre of Excellence in Cambridge (UK).

T: +44 1304 640538
E: charlie.bowen@pfizer.com



Dr Jay Sayed

Jay Sayed, PhD, Principal Device Development Scientist at Pfizer, is a medical engineer with over 10 years of experience in drug-delivery device development. Based at Pfizer's Devices Centre of Excellence in Cambridge, Dr Sayed is focused on the development and manufacture of spring-driven autoinjectors.

T: +44 7590 537503
E: jay.sayed@pfizer.com



Harriet Field

Harriet Field, Laboratory Science Apprentice, has been working at Pfizer since 2021 through the Level 6 Laboratory Scientist apprenticeship programme. Ms Field joined the Device Centre of Excellence in 2024 and previously worked in Chemical R&D. As part of her apprenticeship, Ms Field is studying part-time towards a BSc chemistry qualification, and is now in her final year.

E: harriet.field@pfizer.com

Pfizer Inc

235 East 42nd Street, New York, NY 10017, United States
www.pfizer.co.uk

WHERE CONTENT MEETS INTELLIGENCE

www.ondrugdelivery.com/subscribe



SANNER GROUP

Sanner is a smart, agile and reliable CDMO, providing a full spectrum of services, from the initial design and development phases to mass production of everything from simple, single-use consumables to complex electromechanical devices.

Our agility, scalability and substantial production capacity allows us to cater to a wide range of specific customer needs - whether for small production runs or high-volume manufacturing.

Supported by our state-of-the-art production facilities across Europe, Asia, and North America, Sanner is a strong partner of choice to support our clients to bring medical devices to their local markets efficiently.



28.500 sqm of state-of-the-art production space with available capacities. 11.000 sqm cleanroom space Class 7 & 8



Production footprint on three continents: Europa, North America, and Asia



Highest quality level, GMP compliant and DIN EN ISO 9001, ISO 15378, & ISO 13485 production space



150 engineers for design and development of devices and packaging

**Contact us today
to start your device
design and production!**

sales@sanner-group.com
www.sanner-group.com

**And don't forget
to visit us at:**

PDA Vienna, #150 hall X2
PODD Boston, #26
CPHI Frankfurt, #0J92 hall 8